

KPL SureBlue Reserve™ TMB Microwell Peroxidase Substrate (1-Component)

Catalog No.	Size
5120-0081 (53-00-01)	100 mL
5120-0082 (53-00-02)	400 mL
5120-0083 (53-00-03)	1000 mL

DESCRIPTION

KPL SureBlue Reserve™ TMB Microwell Peroxidase Substrate (1-Component) contains 3,3',5,5' tetramethylbenzidine in a mildly acidic buffer. It develops a deep blue soluble product when reacted with horseradish peroxidase labeled conjugates in microwell plates or tubes. KPL SureBlue Reserve™ provides more sensitivity compared to standard KPL SureBlue TMB substrate. It is suitable for qualitative or quantitative microwell immunoassays and not recommended for membrane or immunohistochemical applications.

FORM/STORAGE/STABILITY

Liquid. Store at 2-8 °C. Do not freeze. Stable for 36 months from date of manufacture when stored at 2-8 °C.

CONTENT

KPL SureBlue Reserve™ TMB Microwell Peroxidase Substrate (1- Component) contains 3,3',5,5'-tetramethylbenzidine in a mildly acidic buffer.

- 5120-0081 (53-00-01) contains 1 x 100 mL
- 5120-0082 (53-00-02) contains 4 x 100 mL
- 5120-0083 (53-00-03) contains 1 x 1000 mL

APPEARANCE/BACKGROUND

The expected appearance of the solution is colorless to very light yellow. SeraCare's QC specification for background of unreacted substrate at 650 nm is ≤ 0.03 Au.

RECOMMENDED HANDLING

KPL SureBlue Reserve™ TMB Microwell Peroxidase Substrate (1- Component) is extremely sensitive to certain handling and storage conditions.

- Avoid exposure to sunlight, air and extreme temperatures, however exposure to indirect laboratory lighting for 8 hours is not harmful. A pale blue color may become apparent after 4 – 5 hours of exposure to indirect laboratory lighting, but this will not affect product performance.
- Never pipette directly from bottle. Pour out needed amount into a plastic reservoir. Do not return excess KPL SureBlue Reserve™ to the primary container.

- Avoid exposure to silica-based materials including borosilicate glass, metals, bacterial contamination or other oxidizing agents.
- Do not refilter as this may destabilize the dye complex and result in the development of background color.
- Redispense into amber Nalgene HDPE and LDPE bottles ONLY. Do not redispense into glass bottles as this may compromise substrate performance.
- HDPE and LDPE bottles should be utilized as they are received directly from the manufacturer; washing of bottles is not recommended.
- Redispense substrate using Platinum Cured Silicone tubing and do not reuse tubing. Avoid redispersing with all other tubings.
- Gravity fill or the use of a peristaltic pump is advised.
- Avoid substrate contact with any metallic surfaces.
- For further details concerning the redispersing of this product, please contact SeraCare's Technical Service Department at 800-638-3167 or 301-948-7755.

USE

Warm to room temperature before use. Solution is ready to use.

Substrate Development

Following incubation with peroxidase labeled conjugate, wash plate thoroughly. Add 100 μ L substrate solution to each well. As the color develops, tap gently to mix. Incubation times will vary depending on your assay.

Over-reaction

Precipitate in the wells after stopping is a sign of over-reaction with KPL SureBlue Reserve™. Precipitate usually takes 10 – 20 minutes to develop so reading the plate immediately after adding the stop solution is recommended. To reduce the intensity of the substrate reaction, further dilution of the primary antibody and/or HRP labeled conjugate is recommended. Dilution of the substrate is not recommended.

To Stop Reaction

For optimal performance, stop reaction by adding equal volume of 1N HCl, 0.6N Sulfuric acid or TMB Stop

KPL SureBlue Reserve™ TMB Microwell Peroxidase Substrate (1-Component)

Catalog No.	Size
5120-0081 (53-00-01)	100 mL
5120-0082 (53-00-02)	400 mL
5120-0083 (53-00-03)	1000 mL

Solution to the microwell plate. This will halt color development and will turn the KPL SureBlue Reserve™ substrate yellow. Additionally, the reaction may be stopped by using KPL BlueSTOP, which retains the blue color.

To Read Reaction

After stopping, read at a wavelength of 450 nm. Stopped reaction should be read within 1 hour.

When to Stop Substrate Reaction

Upon addition of stop solution, absorbance values increase 2 - 3 fold. The point at which the substrate reaction is stopped is often determined by the ELISA reader. The O.D. values of the plate should be monitored and the reaction stopped before positive wells are no longer readable.

ABSORBANCE MEASUREMENTS

Kinetic Assays

The KPL SureBlue Reserve™ substrate produces a blue color upon reaction with peroxidase. Read at a wavelength of 620 - 650 nm.

Endpoint Assays

The addition of 100 µL (or an equal volume) of stop solution to the microwell plate will halt color development and will turn the KPL SureBlue Reserve™ substrate yellow. Read at a wavelength of 450 nm. Stopped reactions should be read within 1 hour.

RECOMMENDED STOP SOLUTIONS

For best results, SeraCare recommends 1N HCl, 0.6 N Sulfuric acid, SeraCare's KPL TMB Stop Solution or SeraCare's KPL TMB BlueSTOP Solution (See Related Products) to stop the reaction.

PRODUCT SAFETY AND HANDLING

See MSDS (Material Safety Data Sheet) for this product.

RELATED PRODUCTS

Cat. No.

RELATED PRODUCTS	Cat. No.
TMB Stop Solution	5150-0020 (50-85-05)
TMB BlueSTOP	5150-0022 (50-85-30)
Wash Solution Concentrate	5150-0008 (50-63-00)
BSA Diluent/Blocking Solution Concentrate	5140-0006 (50-61-00)
Coating Solution Concentrate	5150-0014 (50-84-00)

SureBlue Reserve™ is a trademark of SeraCare, Inc.